


EMERGENCY DEPARTMENT

FOR Securing of Chest Drains


1. Hollister Vertical Tube Attachment device

USED FOR 28 Fr TUBES
(PLUS A HOLDING SUTURE)

Pre-attachment



Post attachment




Consider placement/direction of secure tab to ensure the tab is not pointing in a direction which could cause discomfort.

2. Unomedical Drain Fix

USED FOR 12 Fr SELDINGER TUBES
(PLUS A HOLDING SUTURE)

Pre-attachment



Post attachment

